



# AriaMx HRM Calibration Plate Kit

**Part Number 5190-7702**

## Protocol

**For Research Use Only. Not for use in diagnostic  
procedures.**

**Version B0, October 2021**



**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2014, 2021

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

5991-9555

## Edition

Version B0, October 2021

Printed in USA

Agilent Technologies, Inc.  
5301 Stevens Creek Blvd  
Santa Clara, CA 95051 USA

## Technical Support

Agilent's world-wide Sales and Support Center contact details for your location can be obtained at [www.agilent.com/en/contact-us/page](http://www.agilent.com/en/contact-us/page)

## Warranty

**The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

---

### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

---

## In this Guide...

This protocol provides instructions for using the AriaMx HRM Calibration Plate Kit.

### **1 Before You Begin**

Make sure that you read and understand the information in this chapter before you start an experiment.

### **2 Protocol**

This chapter provides instructions on running an HRM calibration experiment using the AriaMx HRM Calibration Plate.

## What's New in Version B0

- Updated contact information for Agilent Technical Support.

# Content

## 1 Before You Begin

Notice to Purchase 8

Overview 9

Kit components 9

Plate storage and handling 9

## 2 Protocol

Prepare the plate 12

Run the HRM calibration experiment 12

If the HRM calibration experiment fails 13





# 1 Before You Begin

Notice to Purchase 8  
Overview 9  
Kit components 9  
Plate storage and handling 9

Make sure that you read and understand the information in this chapter before you start an experiment.



## Notice to Purchase

This product is provided under an agreement between Biotium, Inc. and Agilent Technologies, Inc. and the manufacture, use, sale or import of this product is the subject to one or more of pending patent applications owned by Biotium and Allelogic. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer, where such research does not include testing, analysis or screening services for any third party in return for compensation on a per test basis. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. Commercial Purposes means any activity by a Party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For information on purchasing a license to this product for purposes other than research, contact Biotium, Inc., Business Development, 3423 Investment Blvd, Suite 8, Hayward, CA 94545 Tel: 510-265-1027, Fax: 510-265-1352.



## Overview

The AriaMx Real-Time PCR System allows you to perform high resolution melt (HRM) analysis. HRM analysis is a technique used for genotyping samples that include a single nucleotide polymorphism (SNP) in the DNA sequence. Applications that may use HRM analysis include species identification, mutation screening, and haplotype characterization.

In order to use HRM analysis on the AriaMx instrument, you must run an HRM calibration experiment using an HRM calibration plate (HCP). The AriaMx HRM Calibration Plate Kit contains a 96-well plate for use in an HRM calibration experiment. The plate is pre-aliquoted with a master mix containing a DNA product and the double-stranded DNA-binding dye EvaGreen.

## Kit components


The AriaMx HRM Calibration Plate Kit contains one 96-well plate. Each of the 96 wells contains 20  $\mu$ L of pre-aliquoted reagents.

## Plate storage and handling

The plate is shipped on dry ice. Store at  $-20^{\circ}\text{C}$  upon receipt.

Wear chemical-resistant powder-free gloves whenever handling the HRM Calibration Plate.

**1 Before You Begin**  
Plate storage and handling



## 2 Protocol

- Prepare the plate 12
- Run the HRM calibration experiment 12
- If the HRM calibration experiment fails 13

This chapter provides instructions on running an HRM calibration experiment using the AriaMx HRM Calibration Plate.



## Prepare the plate

Prepare the plate just before starting an HRM calibration experiment.

- 1 Thaw the AriaMx HRM Calibration Plate at room temperature for 5-15 minutes. *Do not remove the seal.*
- 2 Once thawed, mix the contents of the wells by completing the following actions five times.
  - a Invert the plate so that the liquid moves to the seal on top of each well.
  - b Re-invert the plate to return the liquid to the bottom of each well.
- 3 After completing five mixing cycles, briefly spin the plate in a plate centrifuge at  $500 \times g$  for at least one minute.
- 4 Verify that no bubbles are present at the bottoms of the wells, beneath the liquid.

While it is critical to remove any bubbles from the bottom of the tubes, small bubbles near the liquid meniscus should not affect the test results.

## Run the HRM calibration experiment

You must set up the HRM calibration experiment directly on the instrument (you cannot set up an HRM calibration experiment from the AriaMx software on your PC).

Verify that the FAM/SYBR optical module is installed in the AriaMx instrument before you begin.

- 1 Load the AriaMx HRM Calibration Plate into the AriaMx instrument and close the instrument door.
- 2 On the instrument Home screen, press the HRM Calibration icon.
- 3 On the subsequent screen, press **Open Default Experiment**.

The default HRM calibration experiment opens to the Plate Setup screen. All wells are set to the Unknown well type and the SYBR dye is selected for target detection in all wells.

The EvaGreen dye used in the AriaMx HRM Calibration Plate is detectable with the FAM/SYBR optical module.

**4** Navigate to the Thermal Profile screen and press **Run Experiment**.

A message box opens on the touchscreen prompting you to save the experiment. Press **OK** in the message box to save the experiment to the HCP folder.

**5** Enter a file name for the experiment and press **Save**.


The instrument starts running the experiment.

**6** After the run, a message box opens on the screen notifying you if the calibration passed or failed.

- If it passed, copy the post-run experiment file to your PC to be used for calibration in experiments with an HRM segment.
- If it failed, see “[If the HRM calibration experiment fails](#)”, below.

## If the HRM calibration experiment fails

If the calibration experiment failed the system's quality check, the touchscreen displays a message box at the end of the run notifying you of the failure. You will also see a warning icon (as shown below) if you open the experiment in the AriaMx program on your PC.



Possible causes of a failed experiment include improper handling of the plate and the presence of bubbles in the bottoms of the plate wells. Review “[Plate storage and handling](#)” on page 9 and “[Prepare the plate](#)” on page 12. If problems persist, contact Agilent [Technical Support](#).

[www.agilent.com](http://www.agilent.com)

## In This Book

This guide contains information to run the AriaMx HRM Calibration Plate Kit protocol.

© Agilent Technologies, Inc. 2014, 2021

Version A1, September 2021



5991-9555



**Agilent Technologies**